

HAULWAY™ Scooter & Power Chair Lift

Ordering Checklist

Please complete the following information before quoting or ordering equipment.

Customer Name			Today's Date	
Customer Address		City		State
Zip				
Phone	Cell Phone	Email		

Mobility Device



3-Wheel Scooter



4-Wheel Scooter



Power Chair

Mobility Device Details

Make _____ Model _____ Weight _____ pounds

Wheelbase (if known) _____ inches Dimensions (L x W) _____ inches _____ inches

Hitch

Does your vehicle currently have a hitch? Yes Class _____ No *Product requires a hitch, please contact a local hitch installer*

Hitch Height (Working distance from ground) _____ inches *If distance is greater than 17.5", drop L-tube is required*

Tongue Weight Capacity

_____ + _____ < 300 pounds = Class II
 Weight of Mobility Device Weight of Lift

_____ + _____ < 500 pounds = Class III
 Weight of Mobility Device Weight of Lift

Ordering Equipment

Follow these easy steps to select your perfect setup

STEP 1: Choose the main unit

- Lite
- Standard
- Heavy Duty (Universal Tray)

STEP 2: Select tray configuration

- 3-Wheel (1 center tray)
- 4-Wheel (2 side trays)
- Deluxe (center and side trays)
- Power Chair (2 wide side trays)
- Power Chair Deluxe (center and wide side trays)
- Universal (Only available on HD)

STEP 3: Select optional equipment

- Standard Swingaway
- Heavy Duty Swingaway
- Scooter Cover
- Power Chair Cover